Instructions for Use

Minisart® NML **CE**₀₁₂₃

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Applications

- Minisart® NML with 5 µm hydrophilic cellulose acetate membrane and methacrylate butadiene styrene (MBS) housing are suitable for particle removal and clarification of aqueous solutions.

Intended Purposes

- Minisart® NML, 5 µm cellulose acetate: 17594-----K, 17594-----GUR, customer-specific item number: 17594-----GJR; all versions with MBS housing.
- Intended purpose: Minisart® NML 5 µm are intended to filter low volume aqueous solutions for their clarification to prevent transmission of particulate contaminants during fluid administration into the patient's body in a professional healthcare environment as accessories to syringes or injection systems with Luer connectors.

Intended patient population

- Minisart® filters are intended for use in patients of all age groups, genders, including vulnerable populations. EO gas-sterilized Minisart® filters are not intended to be used for premature babies, stem cells, or embryos.

Clinical benefit

- The clinical benefit of Minisart® NML 5 µm syringe filters is the removal of potential contaminants (i.e., particulate) as part of clarification from aqueous liquids, which are to be administered to the patients in a therapeutic context in a healthcare setting.

Instructions for Use

Minisart® syringe filters may only be used by qualified specialist personnel in healthcare institutions.

1. Filtration during the Preparation of Medication

If possible, draw at least 1 mL air into a syringe (a). Then aspirate the liquid to be filtered from a suitable container into the syringe (see note 1) (b). Open sterile blister units by peeling off the protective backing. Connect the syringe to the Minisart® syringe filter inside the package, while firmly holding the outside of the opened blister packaging. If required, remove excess liquid from the tip of the syringe before attaching it to the Minisart® by Luer lock or Luer slip connection (see note 3) (c). Depress the plunger of the syringe with constant pressure to pass the contents of the syringe through the Minisart® syringe filter into a suitable vial for collection (see warning 6) (d). Afterwards, depress the plunger fully, so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter, to reduce the dead volume. Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).

2. Preparation of Injection Solutions

Draw the liquid to be filtered from a suitable container into the syringe (see note 1) (b). Open sterile blister units by peeling off the protective backing. Connect the syringe to the Minisart® syringe filter inside the package, while firmly holding the outside of the opened blister packaging. If required, remove excess liquid from the tip of the syringe before attaching it to the Minisart® by Luer lock or Luer slip connection (see note 3) (c). Depress the plunger of the syringe until liquid flows out of the Minisart® and ensure that all air is removed before the liquid is injected through the filter into the patient's body (see note 2). Depress the plunger of the syringe with constant pressure to pass the contents of the syringe through the Minisart® syringe filter. Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).

3. Filtration when Aspirating Solutions

Open sterile blister units by peeling off the protective backing. Connect the sterile Minisart® inside to the syringe via the Luer lock or Luer slip connection, while firmly holding the outside of the opened blister package (see note 3). Aspirate the solution to be filtered from a suitable container through the Minisart® into the syringe (see note 1). Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).

Warnings!

- Before use, please check that the Minisart® syringe filter is not damaged and ensure that the unit contains a white membrane to guarantee particle separation according to given pore size. (Note: The filter membrane appears white in color and can be seen inside the plastic housing.) If using sterile Minisart® in blister packaging, please check that the blister packaging is not damaged and use the Minisart® immediately after opening the blister to ensure the sterility of the device. Do not use if you identify any external or internal damage or impurities!
- Minisart® syringe filters are designed for transient (<1 hour) single-use. To help avoid cross contamination, do not use more than once or for more than one sample.
- Do not use Minisart® syringe filters outside of the technical specifications (see Table 1).
- Do not use Minisart® syringe filters that have been sterilized using EO gas for premature babies, stem cells, or embryos.
- Typical filtration volumes of liquids are around 5–100 mL. An accurate prediction of possible filter volumes is not possible due to the different particle load of solutions. A decreasing flow rate at normal effort is clear evidence for approaching maximum filter capacity. To prevent the membrane bursting and guarantee particle separation according to given pore size, do not exceed the maximum operating pressure of 4.5 bar.
- When using syringes with a volume of less than 10 mL, the maximum operating pressure of Minisart® syringe filters of 4.5 bar (65 psi) can quickly be exceeded. Excess pressure could damage the filter, meaning that clarification is no longer guaranteed. We recommend using syringes larger than 10 mL. For syringes smaller than 10 mL, please depress the plunger slowly and if clear resistance is reached (=maximum filter capacity), please do not continue filtration or apply more force, otherwise the filtration pressure will exceed the load limit of the Minisart® syringe filter.
- Minisart® syringe filters can be connected to an active medical device of class Ia or above (for example, a medication pump). Ensure that the pressure in the system is not higher than the maximum operating pressure of the Minisart® syringe filter of 4.5 bar (65 psi). The specific application must be checked by the user in advance.
- Minisart® NML syringe filters are designed for bidirectional use. However, once you have selected a direction of filtration, this direction must be maintained. Do not use the same syringe filter for filtration in both directions!
- Wear suitable protective clothing and protective goggles when using Minisart® syringe filters.
- Do neither use non-sterile Minisart® syringe filters nor those with a pore size of 0.45 µm or 5 µm for sterile filtration. The color coding of the filter makes it possible to confirm which membrane is in the Minisart® syringe filter. A detailed overview of the color coding can be found in the ordering information in Table 2.
- The user must ensure that the filter used is suitable for its application, i.e. the filtered solution is not impaired by the filtration step and only negligible amounts of pharmaceutical active substance are lost. The Minisart® syringe filter is neither suitable for air nor gas filtration.
- Filtration of incompatible media may impair the filter integrity, meaning that it is not possible to ensure particle separation according to the specified pore size. The compatibility of the substances listed in Table 3 has been confirmed either via literature research or tests in the laboratory. However, as the chemical compatibility depends on various factors, the user should check the compatibility of the solution to be filtered before application.
- Certain incompatible drug-drug, drug-solvent and drug-device interactions may lead to particle formation, filter occlusion, and / or low drug dose recovery in the filtrate. The absence of incompatibilities must be validated by the user or organization before using the Minisart® syringe filters.
- Only connect the Minisart® syringe filter using Luer connections. Only this connection type will guarantee a secure connection with other products using the standardized Luer dimensions.
- Store the Minisart® between 5°C and 30°C in a frost-free location and do not subject the Minisart® to direct sunlight or moisture. The medium to be filtered must not exceed the specified maximum temperature (see Table 1). To ensure the integrity and functionality, discard the Minisart® and, if applicable, the filtrate in case the Minisart® was stored or used outside the allowable storage or use conditions.
- If malfunctions occur during use, stop the filtration process immediately and discard the filtrate.

Note

- If you need to spike a closure first before injecting the solution, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile cannula to the outlet of the filter or the syringe and remove the protective cap before injection.
- Hold the syringe with filter and needle pointing upwards. Slowly press the liquid upwards. This assures proper wetting of the membrane filter and prevents air bubbles.
- A tight connection is achieved with half a turn or by plugging into the inlet connection with normal force.
- If you are working with radioactive, toxic, or infectious material, discard the Minisart® syringe filter after use and collect the waste according to the local regulations.
- Minisart® NML, sterile, with a cellulose acetate membrane are sterilized using validated processes with the following parameters: Gamma sterilization: 45 kGy (max.), EO sterilization: Temperature <50°C | relative humidity <85%. The function of the syringe filter was tested by Sartorius using these parameters.
- Do not re-sterilize Minisart® that are supplied sterile in blister packaging.
- All incidents and complaints that occur in relation to the product must be reported to Sartorius or the distributor from which you purchased the product and to the competent authority of the Member State in which the user and / or patient is established.

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3035372 /MAN / 000 / 00 / Applicable Site: / PDF Creation Date: 2025-04-03

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Instructions for Use

Minisart® NML CE0123

Sartorius Single-use Syringe Filters for the Filtration of Aqueous Solutions


Applications
• Minisart® NML with 5 µm hydrophilic cellulose acetate membrane and methacrylate butadiene styrene (MBS) housing are suitable for particle removal and clarification of aqueous solutions.

Intended Purposes
• Minisart® NML, 5 µm cellulose acetate: 17594-----K, 17594-----GJR, all versions with MBS housing, customer-specific item number: 17594-----GJR; all versions with MBS housing.
Intended purpose: Minisart® NML 5 µm are intended to filter low volume aqueous solutions for their clarification to prevent transmission of particulate contaminants during fluid administration into the patient's body in a professional healthcare environment as accessories to syringes or injection systems with Luer connectors.

Intended patient population
• Minisart® filters are intended for use in patients of all age groups, genders, including vulnerable populations.
EO gas-sterilized Minisart® filters are not intended to be used for premature babies, stem cells, or embryos.

Clinical benefit
• The clinical benefit of Minisart® NML 5 µm syringe filters is the removal of potential contaminants (i.e., particulate) as part of clarification from aqueous liquids, which are to be administered to the patients in a therapeutic context in a healthcare setting.

Instructions for Use
Minisart® syringe filters may only be used by qualified specialist personnel in healthcare institutions.

- Filtration during the Preparation of Medication**
If possible, draw at least 1 mL air into a syringe (a). Then aspirate the liquid to be filtered from a suitable container into the syringe (see note 1) (b). Open sterile blister units by peeling off the protective backing. Connect the syringe to the Minisart® syringe filter inside the package, while firmly holding the outside of the opened blister packaging. If required, remove excess liquid from the tip of the syringe before attaching it to the Minisart® by Luer lock or Luer slip connection (see note 3) (c). Depress the plunger of the syringe with constant pressure to pass the contents of the syringe through the Minisart® syringe filter into a suitable vessel for collection (see warning d) (d). Afterwards, depress the plunger fully, so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter, to reduce the dead volume. Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).
- Preparation of Injection Solutions**
Draw the liquid to be filtered from a suitable container into the syringe (see note 1) (b). Open sterile blister units by peeling off the protective backing. Connect the syringe to the Minisart® syringe filter inside the package, while firmly holding the outside of the opened blister packaging. If required, remove excess liquid from the tip of the syringe before attaching it to the Minisart® by Luer lock or Luer slip connection (see note 3) (c). Depress the plunger of the syringe until liquid flows out of the Minisart® and ensure that all air is removed before the liquid is injected through the filter into the patient's body (see note 2). Depress the plunger of the syringe with constant pressure to pass the contents of the syringe through the Minisart® syringe filter. Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).
- Filtration when Aspirating Solutions**
Open sterile blister units by peeling off the protective backing. Connect the sterile Minisart® inside to the syringe via the Luer lock or Luer slip connection, while firmly holding the outside of the opened blister package (see note 3). Aspirate the solution to be filtered from a suitable container through the Minisart® into the syringe (see note 1). Remove the Minisart® and dispose of it in accordance with the applicable regulations (see note 4).

- Warnings!**
- Before use, please check that the Minisart® syringe filter is not damaged and ensure that the unit contains a filter membrane to guarantee particle separation according to given pore size. (Note: The filter membrane appears white in color and can be seen inside the plastic housing.) If using sterile Minisart® in blister packaging, please check that the blister packaging is not damaged and use the Minisart® immediately after opening the blister to ensure the sterility of the device. Do not use if you identify any external or internal damage or impurities!
 - Minisart® syringe filters are designed for transient (< 1 hour) single-use. To help avoid cross contamination, do not use more than once or for more than one sample.
 - Do not use Minisart® syringe filters outside of the technical specifications (see Table 1).
 - Do not use Minisart® syringe filters that have been sterilized using EO gas for premature babies, stem cells, or embryos.
 - Typical filtration volumes of liquids are around 5-100 mL. An accurate prediction of possible filter volumes is not possible due to the different particle load of solutions. A decreasing flow rate at normal effort is clear evidence for approaching maximum filter capacity. To prevent the membrane bursting and guarantee particle separation according to given pore size, do not exceed the maximum operating pressure of 4.5 bar.
 - When using syringes with a volume of less than 10 mL, the maximum operating pressure of Minisart® syringe filters of 4.5 bar (65 psi) can quickly be exceeded. Excess pressure could damage the filter, meaning that clarification is no longer guaranteed. We recommend using syringes larger than 10 mL. For syringes smaller than 10 mL, please depress the plunger slowly and if clear resistance is reached (=maximum filter capacity), please do not continue filtration or apply more force, otherwise the filtration pressure will exceed the load limit of the Minisart® syringe filter.
 - Minisart® syringe filters can be connected to an active medical device of class IIa or above (for example, a medication pump). Ensure that the pressure in the system is not higher than the maximum operating pressure of the Minisart® syringe filter of 4.5 bar (65 psi). The specific application must be checked by the user in advance.
 - Minisart® NML syringe filters are designed for bidirectional use. However, once you have selected a direction of filtration, this direction must be maintained. Do not use the same syringe filter for filtration in both directions!
 - Wear suitable protective clothing and protective goggles when using Minisart® syringe filters.
 - Do neither use non-sterile Minisart® syringe filters nor those with a pore size of 0.45 µm or 5 µm for sterile filtration. The color coding of the filter makes it possible to confirm which membrane is in the Minisart® syringe filter. A detailed overview of the color coding can be found in the ordering information in Table 2.
 - The user must ensure that the filter used is suitable for its application, i.e. the filtered solution is not impaired by the filtration step and only negligible amounts of pharmaceutical active substance are lost. The Minisart® syringe filter is neither suitable for air nor gas filtration.
 - Filtration of incompatible media may impair the filter integrity, meaning that it is not possible to ensure particle separation according to the specified pore size. The compatibility of the substances listed in Table 3 has been confirmed either via literature research or tests in the laboratory. However, as the chemical compatibility depends on various factors, the user should check the compatibility of the solution to be filtered before application.
 - Certain incompatible drug-drug, drug-solvent and drug-device interactions may lead to particle formation, filter occlusion, and | or low drug dose recovery in the filtrate. The absence of incompatibilities must be validated by the user or organization before using the Minisart® syringe filters.
 - Only connect the Minisart® syringe filter using Luer connections. Only this connection type will guarantee a secure connection with other products using the standardized Luer dimensions.
 - Store the Minisart® between 5°C and 30°C in a frost-free location and do not subject the Minisart® to direct sunlight or moisture. The medium to be filtered must not exceed the specified maximum temperature (see Table 1). To ensure the integrity and functionality, discard the Minisart® and, if applicable, the filtrate in case the Minisart® was stored or used outside the allowable storage or use conditions.
 - If malfunctions occur during use, stop the filtration process immediately and discard the filtrate.

- Note**
- If you need to spike a closure first before injecting the solution, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile cannula to the outlet of the filter or the syringe and remove the protective cap before injection.
 - Hold the syringe with filter and needle pointing upwards. Slowly press the liquid upwards. This assures proper wetting of the membrane filter and prevents air bubbles.
 - A tight connection is achieved with half a turn or by plugging into the inlet connection with normal force.
 - If you are working with radioactive, toxic, or infectious material, discard the Minisart® syringe filter after use and collect the waste according to the valid regulations.
 - Minisart® NML, sterile, with a cellulose acetate membrane are sterilized using validated processes with the following parameters: Gamma sterilization: 45 kGy (max.), EO sterilization: Temperature <50°C | relative humidity <85%. The function of the syringe filter was tested by Sartorius using these parameters.
 - Do not re-sterilize Minisart® that are supplied sterile in blister packaging.
 - All incidents and complaints that occur in relation to the product must be reported to Sartorius or the distributor from which you purchased the product and to the competent authority of the Member State in which the user and | or patient is established.

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© 2025 Sartorius Stedim Biotech GmbH
Sartorius Australia Pty. Ltd.
Sartorius Stedim Biotech GmbH